

This is an EXAMPLE COMPRESSION test method. This is a PROMPTED TEST where you are "prompted" step-by-step.
This example method is "Read-Only".

Operator ID	George
Laboratory Name	OlinMatSci
Company	EEF
Temperature (deg C)	18.00
Humidity (%)	50.00
Rate 1	10.00000 mm/min
Number of specimens	3



Graph 1



	Specimen label	Maximum Load (kgf)	Yield Strength (Offset 0.2 %) (MPa)	Compressive Strength (MPa)	Comment
1	1	0.73	0.01	0.02	
2	2	0.42	0.00	0.01	
3	3	0.56	0.01	0.01	
4	4	0.35	0.00	0.01	
5	5	0.51	0.01	0.01	interesting side diagonal breaks
6	6	0.40		0.01	
7	7	0.44		0.01	
8	8	0.72		0.02	
Mean		0.52	0.01	0.01	
Standard Deviation		0.14200	0.00170	0.00354	
Minimum		0.35	0.00	0.01	
Maximum		0.73	0.01	0.02	
Range		0.38	0.00	0.01	